

List with tested cleaning products:

> Always follow the prescribed concentration

> Always rinse with water and neutral soap after cleaning

Using cleaning agents can cause discolouration, but there is no deterioration of the coating layer.

> Products pre-marked with a "V" are listed on the VAH list of disinfectants

Type	Surface disinfection		Wound disinfection (100%)
Alcohol	v B40	100%	Kodan® Tinktur Forte farblos
	v Bacillo® 30 Foam	100%	
	v Bacillo® 30 Tissues	Wipes	
	Dax Ytdesinfektion Plus	100%	
	Ethanol	75%	
	Isopropanol	45%	
	v Meliseptol®	100%	
	v Mikrozyd® AF Liquid	100%	
	v Mikrozyd® AF Wipes	Wipes	
	v Mikrozyd® universal wipes premium	Wipes	
	v RUCK Sprühdeseinfektion	100%	
	Sagrotan® Hygiene-spray	100%	
	Sprint 200	2%	
Ammonia	Volvone	40%	
Chlorine	Actichlor Plus	5 000 ppm	Hansamed spray
	Chlor-Clean tablets	5 000 ppm	Hibicet
	Germatol	5%	Octenisept®
	Haz-Tabs	5 000 ppm	
	Sodium hypochlorite	5 000 ppm	
Hydrogenperoxide	Hydrogenperoxide	6%	
	v Incidin® OxyFoam	100%	
	v Incidin® OxyFoam S	100%	
	v Incidin® OxyWipe S	Wipes	
	Wecoline Clean 'n Easy	100%	
	Laquick Pro	100%	
Oxygen	v Anios Oxy-Floor	2%	
	v Dismozon® pur	6%	
	Virkon Tablets	1%	
Peracetic Acid	Clinell Sporidical Wipes	Wipes	
	Terralin® PAA (pH=5,2)	6%	
	v Ultrasol® active	5%	
Quaternary	v Biguanid Fläche N	4%	
	Incidin Extra N	5%	
	RUCK Maxitücher Alkoholfrei	Wipes	
	RUCK Schaumdesinfektion	100%	
	Suma Bac D10	1%	
	Surfanios Premium	3%	
	v SurfaSafe Premium	100%	
	v Terralin® protect	4%	
	Umonium medial spray	100%	
Quaternary & Alcohol	v B 30	100%	
	v B 45	100%	
	v Incidin® Foam	100%	
	Incidin® Rapid	2%	
	Incidin® Pro	5%	
	v RUCK Flachendesinfektionstucher, alkoholisch	Wipes	
	v RUCK Flachendesinfektionstucher, alkoholfrei	Wipes	
Other	Alkaline		
	Alkylamine	v P3-topax® 66	5%
	Aliphatic Hydrocarbons	v Incidin® Plus	4%
	Amino alcohol	RUCK Kunststoffreiniger	100%
	Chloroxylenol	Sanit P20	100%
	Formaldehyde	Dettol	5%
	Glycol	Formaldehyde	30%
	Sodium hydroxide	Spring	100%
	Topmax 421	5%	

Type	Cleaning		
Surfactants	Reocid (pH=5,6) Universal wipes Dasty Classic degreaser (pH=11.2) Biosafe D	5% Wipes 100% 10%	
	Cleaner 500	100%	

AD2 Textiles | Loodstraat 34 | 2718 RW ZOETERMEER | The Netherlands